


How to attach the Click Hock Shield Ultra

VERY IMPORTANT: Make sure the leg is straight under and the hoof is flat on the ground.

1. Make sure the side with tag is on the outside of leg.
2. Pull front flap around to outside of leg. Pull it only so it will connect with the Hock Shield. Have it as loose as it can be without creating wrinkles under the strap.
3. Lay Hock Shield on top of flap attaching inside Velcro pieces.
4. Bring the elastic strap around the inside of leg to the outside and attach to upper Velcro for a “snug” fit.
5. Lock top Velcro flap down.



Step 1



Step 2



Step 3



Step 4



Step 5