


Med Cup Colour Code Guide



- The med cup is colour coded with a Date of Manufacture Label on the connector
- The label can be coloured Grey/Clear, Green or Blue
- The Grey/Clear label med cup is the standard cup that will nebulise most medications
- A Green labelled cup is special for shorter delivery time with some medications
- A Blue labelled cup is special for use with the natural therapy Nebulin®

Important Information:

The med cup & drug combination choice can have a big impact on the delivery time and aerosol droplet size.

- A GREY/CLEAR labelled cup used with high viscosity oily type drugs can lead to slow delivery or failure to nebulise.
- A GREEN labelled cup used with low viscosity watery type drugs can lead to very fast delivery and the generation of big aerosol droplets leading to pools of medication wasted in the bottom of the mask chamber.
- A BLUE labelled cup used with low viscosity watery type drugs will nebulise but the delivery time will be slow. It is not suitable for high viscosity oily type
- If the drug you wish to nebulise is not listed on the drug table then it is best to test a sample of the drug with your cup first to check if it nebulises ok and time before delivering the treatment.
- If you are unsure about what med cup to use it is best to contact your local distributor for advice.
- Check out our website at www.nortev.com/flexineb/med-cup-choices for any updates to this information.

It is recommended to always consult your veterinarian for advice regarding diagnosis and the treatment protocol for your horse.

Antibiotics / Antifungal					
Drug	Trade Name	Drug Format	Dilute with Saline	Standard Med Cup Colour Code *	Special Med Cup Colour Code For Shorter Delivery Time *
Cefquinome	Cobactan®	Cobactan® 4.5% Injectable Solution	No	GREY	GREEN
Ceftiofur Sodium	Excenel®	Excenel® 4G Injectable Solution	No	GREY	GREEN
Gentamicin	Gentaject	Gentaject 10% Injectable Solution	Yes 1:1	GREY	GREEN
Amikacin	Amikavet	Amikavet Injectable Solution	Yes 1:1	GREY	GREEN
Marbofloxacin	Marbocyl	Injectable Solution	Yes 1:1	GREY	GREEN
Nystatin	R-Nystatine	Oral Solution	Yes 1:1	GREY	GREEN
Bronchodilators					
Drug	Trade Name	Drug Format	Dilute with Saline	Standard Med Cup Colour Code *	
Ipratropium	Atrovent®	Atrovent® 250 UDVS Nebuliser solution	Yes 1:1	GREY	
Salbutamol (Albuterol)	Ventolin™	Ventolin™ 2.5 Nebuliser Solution	Yes 1:1	GREY	
Clenbuterol	Ventipulmin®	Ventipulmin® 30 Injectable Solution	Yes 1:1	GREY	
Corticosteroids					
Drug	Trade Name	Drug Format	Dilute with Saline	Standard Med Cup Colour Code *	
Dexamethasone	Dexameth	Dexameth 2mg/ml Injectable Solution	Yes 1:1	GREY	
Flucicasona Propionate	Flixotide™	Flixotide™ 2mg/2ml Nebuliser Suspension	Yes 1:1	GREY	
Budesonide	Pulmicort®	Pulmicort® 0.5 Nebuliser Solution	Yes 1:1	GREY	
Mucolytics					
Drug	Trade Name	Drug Format	Dilute with Saline	Standard Med Cup Colour Code *	Special Med Cup Colour Code For Shorter Delivery Time *
Acetylcysteine	Parvolex®	Parvolex® 2g Solution for infusion	Yes 1:1	GREY	GREEN
0.9% Saline Solution	0.9% Saline Solution	B.Braun NaCl 0.9%	No	GREY	None
Bromhexine	Bisolvon	Oral Solution	Yes 1:1	GREY	None
Emser Salt	Emser Salt	Nebuliser Salt Solution	No	GREY	None
Other Therapies					
Active Ingredients	Trade Name	Format	Dilute with 0.9% Saline	Standard Med Cup Colour Code *	Special Med Cup Colour Code For Shorter Delivery Time *
Natural hydro-soluble extract of curcuma longa	Nebulin®	Powder + NaCl	No		None
Cholated silver and distilled water	SilvaPlex™ & EquiSilver™	Oral Solution	No	GREY	None
Thymus proteins, zinc gluconate, silver solution	BIOPRO-PLUS™	Oral Solution	Yes 1:1	GREY	None
Hyaluronic acid, eucalyptus essential oil in isotonic solution	Hyalusol®	Nebuliser Solution	No	GREY	None
Natural solutions including honey and bryonia cretica	Laryx/Apts	Oral Solution	Yes 1:1	GREY	None
Mixture of essential oils	Inhal-14	Nebuliser Solution	Yes 1:1	GREY	GREEN
Grindella, water and essential oils	Balsamic Air	Oral Solution	Yes 1:1	GREY	GREEN
Urtica dioica, vitis vinefera, avena sativa	Epistop-L	Nebuliser Solution	Yes 1:1	GREEN	None
Mixture of essential oils	StopBlood	Oral Solution	Yes 1:1	GREEN	None

* Read the reverse side of this Drug Table for important information and advice on Med Cup & Drug combination choices

Flexineb Single Medication Cup by Flexineb.

This spare medication cup is used for dispensing a singular medication and should be washed well after each use.

Product Details:

Performance of the medication cup is dependent on the drugs and therapies nebulized and how well it is maintained between uses.

Wash well after each use

Do not immerse the cup in water or liquid

Read product labeling for full cleaning instructions.

Grey/Clear med cup is the standard cup that will nebulise most medications

Grey/Clear cup used with high viscosity oily type drugs can lead to slow delivery or failure to nebulise

Green med cup is special for shorter delivery time with some medications

Green cup used with low viscosity watery type drugs can lead to very fast delivery and the generation of big aerosol droplets leading to pools of medication wasted in the bottom of the mask chamber

Blue med cup is special for use with the natural therapy Nebulin

Blue cup used with low viscosity watery type drugs will nebulise but the delivery time will be slow -

Blue med cup is not suitable for high viscosity oily type drugs

Maintenance / Cleaning Medication Cup: Cleaning: It is very important to wash out the Med Cup after each use to avoid contamination or damage of the metallic vibrating mesh. Protect the electrical connector from getting wet when cleaning by using the protective cap provided. To clean out the Med Cup open up the soft blue Cap and place 2-3 drops of mild dishwashing liquid (do not use dishwasher detergent) into the Cup and fill with distilled water. Close the soft blue Cap and shake the Cup vigorously for about 15 seconds and squeeze the soft Cap to push water through the mesh whilst shaking. Open the Cap and rinse out from both sides with distilled water. Shake off excess water and let air dry before re-use.

Disinfection: To disinfect the outside surfaces of the Med Cup use anti-bacterial wipes.

CAUTION:

Do not immerse the Med Cup in water or cleaning chemical solutions!

Do not nebulise or put chemical cleaning solutions inside the Cup as they may permanently damage the metallic mesh or reduce its life!