

Sizing

Use the size guide to ensure you have the closest fit for your horse or pony. Note - If your horse has specialist or remedial shoes, or large shoes such as Eggbars or Heartbars, you may need to go up a size to accommodate these. We recommend working with your farrier to ensure you get the most applicable size for your horse.

If it is your intention to use on both front, or hind feet, and your horse is particularly 'active' or a bit close in his action, the risk of standing on the boot with the opposite foot is higher due to the decreased distance between the two feet. In this case you can use an overreach boot over the top of the Medical Hoof Boot to protect it.

Size	Width of Hoof (mm)	Width of Hoof (in) approximate
0	76-85	3-3.3
1	86-94	3.4-3.7
2	95-103	3.75-4
3	104 -112	4.1-4.4
4	113-121	4.5-4.75
5	122-130	4.8-5.1
6	131-139	5.2-5.5
7	140-148	5.5-5.8
8	149-157	5.8-6.2
9	158-166	6.2-6.5

