



## Boot Size Guide

|    |               |      |
|----|---------------|------|
| A. | Full Length   | 14"  |
| B. | Bottom Length | 7"   |
| C. | Top Length    | 7"   |
| D. | Top Width     | 8"   |
| E. | Mid Width     | 6"   |
| F. | Bottom Width  | 5.5" |

## How to attach the Click Knee Shield

1. Attach bottom strap just above the fetlock joint. Pull strap snug enough so it doesn't slip over the joint. **DO NOT PULL STRAP TIGHT!**
2. Attach the remaining Velcro tabs loosely.
3. The Knee Shield should fit loosely around the knee. You should be able to slip your fingers between the shield and the knee joint.



**Step 1**



**Step 2**



**Step 3**