

Easyboot REMEDY

Sizing Information

Measuring and sizing for the Easyboot Remedy is different than the other hoof boots in the EasyCare boot range. Typically, boots need a snug fit to prevent twisting and turning. The Easyboot Remedy is different. With the Easyboot Remedy the fit should be loose to allow room around the hoof and lower leg. The extra room allows for the liquid and/or ice to circulate.

The Easyboot Remedy is a loose-fitting boot designed to hold liquid while a horse is tied and monitored.

Easyboot Remedy Size Chart		
	Hoof Width Inches	Hoof Width Millimeters
Easyboot Remedy Size Medium	Hoof width up to 4 3/4 Inches	Hoof width up to 120.5mm
Easyboot Remedy Size Large	Hoof width up to 6 1/2 Inches	Hoof width up to 160mm



EasyCare Inc.
The Ultimate in Hoof Protection