

Novaguard Sizing Charts

Use these 2 charts to size Regular Snout dogs

If you are between 2 sizes, then pick the smaller size.


Steps: 1) Head Measurement 2) Then Canine Weight

1) HEAD MEASUREMENT

(Reference Inches or Centimeters)

< 11	11 - 13	13 - 15	15 - 18	18 - 21	21 - 24	24 +	Centimeters
< 4.3	4.3 - 5	5 - 6	6 - 7	7 - 8.3	8.3 - 9.4	9.4 +	Inches

NOVAGUARD™ SIZE	Mini						
	Toy						
	Small						
	Small Medium						
	Medium						
	Large						
	X-Large						



Measure from **Nose Tip** to **Mid Point** behind ears.

2) CANINE WEIGHT

(Reference Weight in Kilograms and Pounds)

< 2	2-6	5-12	10-25	23-35	33-45	45 +	Kilograms
< 4	4-13	11-26	22-55	50-77	72-99	99+	Pounds

NOVAGUARD™ SIZE	Mini						
	Toy						
	Small						
	Small Medium						
	Medium						
	Large						
	X-Large						

**** Our sizing chart is a guide only ****

**** It does not constitute a guarantee as to the fit of the Optimizor/Novaguard ****