

For: Aid in the temporary relief of minor stiffness and soreness caused by overexertion. Ice-O-Gel cools the leg and tightens at the same time.

Directions: Apply freely to leg first rubbing against hair and then apply more with the hair. Wrap with paper towel and bandage, leaving the treatment on the leg 48 to 72 hours. Repeat if necessary. Leave Ice-O-Gel on leg while working or walking. Remove with warm water. To freeze horse before racing, apply Ice-O-Gel 24 to 36 hours prior. Before post time, reapply Ice-O-Gel using wet, cold water bandages. Place leg in ice tub or ice boots.

Contents: Alcohol.....23.9%
Menthol.....1.0%
Camphor1.2%

Caution: Keep away from children. For veterinary use only. Do not use on food animals.
Store at room temperature. Stir well before using if ingredients are separated.



ICE-O-GEL®

TIGHTENER/FREEZE

Causes little or no scurfing • Cools, draws, and tightens

Net Weight 16 U.S. oz. (.473 Liter)

NDC #049343-9910
DIN 00455431

Hawthorne Products
16828 N. State Road 167 N.
Dunkirk, IN 47336 • 1-800-548-5658
www.hawthorne-products.com

MADE IN USA

