

**Xzantrol
Equine Nasal Flush
& Decongestant**

INDICATIONS:

For the relief of sinus and nasal congestion in the equine. For veterinary use only.

DIRECTIONS:

Shake well before using. Using dosage syringe with catheter tube attached insert tube into bottle of Xzantrol and withdraw 10 cc's of solution into tube and syringe. Slide tube up nostril until only the syringe is visible. Forcefully depress plunger and withdraw tube from nostril. Store at room temperature 59-86 F.

NOTE:

In horses with sever mucous congestion nasal passages may drip for several days after application.

RECOMMENDED APPLICATIONS:

For horses with mild to moderate nasal congestion Xzantrol only needs to be used approximately once every three weeks. For moderate to severe cases Xzantrol can be administered weekly if needed.

CAUTION:

Avoid contact with eyes. Keep this and all medication out of reach of children. Discontinue use and consult veterinarian if condition worsens or fails to improve.

INGREDIENTS:

Aceatic Acid, Potassium Iodide, Saline Solution, Eucalyptus, Camphor Spirits, Natural & Artificial Flavors.

IMPORTANT:

Do not use a veterinary guttural pouch rod as it goes up too high. Xzantrol works best at length of tube supplied. Rinse catheter tube thoroughly with water after withdrawin from each nostril. Failure to rinse tube after withdrawing from each nostril before inserting tube back into Xzantrol solution will result in small mucous particles floating within the solution. Upon completion of administering Xzantrol to both nostrils rinse inside and outside catheter tube thoroughly with hot water or isopropyl alcohol.