

"Silver promotes rapid healing"

EquiSilver™

Veterinary Antimicrobial Silver

Respiratory Solution

100% NATURAL AND NON-TOXIC
16 fl. oz (473 ml)

For use in equine respiratory therapy. Chelated silver inhibits bacterial and viral infections in the sinus cavity or lungs. Equisilver has shown to inhibit low-grade microbial infections in the respiratory tract.

Protocol: Add 4oz of Equisilver to 12oz of distilled water per treatment, continue as needed. Flexineb protocol: Add 25cc of undiluted Equisilver, nebulize 2x daily for 3 days, continue 25cc daily as needed.

*Consult your Veterinarian if problems persist.
Active ingredients: Exclusive patented formulation - Chelated silver and distilled water.

502-387-4202

www.equisilver.com

